

Mini-Pearls™

Compounding System

The **Mini-Pearls™ Compounding System** is excellent for making a high volume of small dose troches or suppositories quickly. It has 228 small spherical cavities with a diameter of 6.3 mm and a volume of 0.13 mL ea. Average Mini-Pearl weight is 0.2 g of Techna 20™ SF (Sweetened/Sugar-Free) Troche Base (SKU: B0070). Durable, flexible, non-stick white silicone food grade top mold that is oven, microwave, freezer, and dishwasher safe. Temp resistant range: -60 to 230°C (-76 to 446°F). Lightweight yet rigid high-density black polymer set plate will not scratch countertops and has no sharp edges: temp resistant range -50°C to 90°C (-58°F to 194°F) and includes four stainless steel alignment pegs. System patent pending.



System Includes:

Item	SKU	Qty
Step-by-step use and cleaning instructions		1
Mini-Pearls™ Silicone Top Mold	E0591-00	1
Mini-Pearls™ Black Set Plate	E0589-00	1
Roller		1
Scraper Spatula, Large (5")	E0569-00	1
Oral Syringe (10 mL)		1
Mini Heating/Sealing Hot Air Tool	E0570-00	1
Multi-Purpose Laboratory Tote with Storage Cover		1
Techna 20™ SF Troche Base (Sweetened)	B0070-08	1 Sample

Recommended Cleaning Method:

Clean by hand with **Relianox™ Critical Cleaning Liquid Detergent** (SKU: P3023).